PL - Stabilizator barku

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Shoulder Stabilizer

USAGE INSTRUCTIONS:

Pass the strap through the closure loop and fasten the Velcro before placing it on the shoulder. Then wrap the cross strap under the opposite arm and secure the Velcro in the center of the chest. Adjust all straps to the desired position, ensuring they are not too tight to avoid discomfort. Next, adjust the pressure pad on the shoulder to fit precisely with the anatomical curves. If possible, ensure that the area supported by the stabilizer is adequately elevated and relieved to minimize the risk of further injury. It is recommended to check the shoulder every 15 minutes. If redness, irritation, or discomfort occurs, remove the stabilizer immediately. If shoulder pain persists, seek a proper diagnosis from a doctor to determine further action.

WARNING:

Do not use the stabilizer while sleeping, near heat sources, or during exercise. Do not apply additional pressure to the stabilizer. Ensure proper air circulation around the stabilizer. If skin discoloration occurs, discontinue use immediately. Not suitable for children, elderly

individuals, or those unable to operate the stabilizer independently. Do not use on irritated or damaged skin. Avoid using the product in combination with creams applied to the skin or in cases of rheumatoid arthritis, poor circulation, nerve damage, sensory disorders, skin conditions, vascular disorders (e.g., Raynaud's disease), diabetes, or pregnancy. For external use only. Do not use the stabilizer if allergic to latex. The product contains natural latex rubber, which may cause allergic skin reactions.

IMPORTANT: To facilitate blood circulation, remove the stabilizer for at least 10 minutes three times a day. Plastic packaging can pose a suffocation hazard—keep it out of reach of children and infants.

WASHING INSTRUCTIONS:

Hand wash only. Do not bleach. Air dry while hanging. Do not iron. Do not dry clean. **MATERIAL:** 75% Neoprene, 25% Polyester

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