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EN - Heel Cushions for Shoes

Heel cushions are intended for use in closed-heel footwear.

Instructions for use: Fit the cushion to the shoe, trim if necessary, remove the protective film, and stick it to the back of the shoe, pressing firmly. Wait a few minutes for the adhesive to take effect.

Warnings: Plastic bags can be dangerous. To avoid the risk of suffocation, keep out of reach of children. If you have foot conditions or pain, consult a doctor or orthopedic specialist before using the cushions. Do not use if the skin is damaged. Avoid using heel cushions if you have open wounds, burns, or irritation around the heel to prevent worsening of the skin condition. Heel cushions are intended for use in closed-heel footwear. Avoid using them in sandals or other types of open footwear, as they may shift and cause discomfort. Do not use if allergic to materials to avoid allergic reactions. Using heel cushions in shoes that are too tight or too loose may cause increased friction, discomfort, or even foot injuries. If the cushions are cracked, torn, or have damaged adhesive surfaces, do not use them as they may shift in the shoe and cause injury. Do not store in inappropriate conditions. Do not leave the cushions in high temperatures (e.g., in the sun, in a car) or in places with high humidity, as this may reduce their effectiveness. Heel cushions are a personal-use product. Do not share them with others to avoid the risk of transmitting bacteria and fungi.

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