Producent/Manufacturer/Hersteller/Valmistaja/Tootja/Ražotājs/Gamintojas/Výrobce/Výrobca/Gyártó/Producător/Производитель
VITOLOG. PL sp. z o.o.
ul. Międzyleska 6
50-514 Wrocław
Nip:8992838703,
Kontakt: +48 71 727 62 28
kontakt@vitolog.pl

Wyprodukowano w Chinach/Made in China/Hergestellt in China/Valmistettu Kiinassa/Valmistatud Hiinas/Ražots Ķīnā/Pagaminta Kinijoje/Vyrobeno v Číně/Vyrobené v Číne/Kínában készült/Fabricat în China/Произведено в Китае

PL:

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User Manual:

Preparation:

Wash and dry the skin in the area where the cups will be applied.
Apply oil or cream to the skin to facilitate cup movement.

Cup Application:

Choose the appropriate cup size depending on the body area.
Connect the pump to the cup using an extension (optional).
Place the cup on the skin and slowly pump until you feel slight resistance and suction.
Leave the cup on the skin for 5-15 minutes depending on individual sensations and therapeutic effects.

Removing the Cup:

To remove the cup, gently press the side and tilt it to allow air in.
If necessary, gently massage the application area.

Aftercare:

After the massage, apply a soothing balm or cream.
Avoid intense physical activity for 2-3 hours after the treatment.

Cleaning and Maintenance:

After each use, wash the cups with warm soapy water and dry thoroughly.
Store the cups in a dry and clean place.

Warnings and Precautions:

Do not apply cups to damaged skin, wounds, bruises, or areas with varicose veins. Avoid using them on the face and near the heart. Do not leave cups on the skin for more than 15 minutes. If you experience discomfort or severe pain, stop the procedure immediately. Do not use on children under 7 years old or pregnant women without consulting a doctor. Store the set out of reach of children. Plastic bags can be dangerous – to avoid suffocation risk, keep them out of reach of children and infants.

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